

Crazy HorsEPower 400R

CRAZY HORSE

CRAZY HORSE **EPOWER**



HMV.Nr.: 18.99.04.00

Version

Electric pulling aid for the wheelchair

- 400 Watt powerful freewheel brushless motor 35 Nm
- Driving function forward and backward
- Speed 6 km/h
- 30 km/h or more for use on private land
- Cruise control

Scope of application



For activities outside the home:

- rolls easily over any surface
- on cobblestones
- sidewalks fall off to one side
- gravel roads, dirt roads, forest roads
- light sleepers and stages
- sandy beach
- shingle beach
- meadows
- city breaks and much more



Steering head

The unique steering head design delivers razor sharp handling and high steering stability and conceded during cornering and straight-line stability, additionally reinforce two springs this effect!



Light Version

400R with hole rims 36 volts - 20 inch rear wheels on 80 mm rims and 80 mm hubs.



Special features



- Frame protection made of plastic
- Hand-parking brake
- 48 volt battery
- Light with horn
- Reverse gear
- Cruise control
- Rotary throttle

