

Crazy HorsEPower 500R

CRAZY HORSE



HMV.Nr.: 18.99.04.00

Version

Electric pulling aid for the wheelchair

- 1100 Watt powerful freewheel brushless motor 65 Nm
- Motor development approx 1000 Watt
- Driving function forward and backward
- Speed 6 km/h
- 15 km/h require
- 30 km/h or more for use on private land
- Cruise control

Special features



- Frame protection made of plastic
- Hand-parking brake
- 48 volt battery
- Light with horn
- Reverse gear
- Cruise control
- Rotary throttle



Scope of application



For activities outside the home:

- rolls easily over any surface
- on cobblestones
- sidewalks fall off to one side
- gravel roads, dirt roads, forest roads
- light sleepers and stages
- sandy beach
- shingle beach
- meadows
- city breaks and much more

