

# Crazy Horse Active Power

CRAZY HORSE

**AKTIV-POWER**



## Version

### Sports Equipment for Wheelchair

- Weight: 7 kg with assembly stand
- Wheel size: 20 inches
- Tyres from 20 x 2.00
- Parking brake may be operated with a single hand
- Disc brake, mechanical, diameter: 160 mm
- Positive steering through special O-ring
- Stand
- Bike bell

## Scope of application



For activities outside the home:

- They roll easily over any surface
- On cobblestones
- Walkways fall off to one side
- Gravel paths, gravel roads, forest roads
- Light-thresholds and steps
- Sand beach (you can sometimes be pushed without immediately falling out of the wheelchair)
- Gravel beach
- Meadows
- Vacation

## Over hill and dale



In our dimension sheet you can assemble your Crazy Horse Active - Power according to your wishes in many variants.

Special orders please inquire!

## Special features



Extremely resilient aluminium frame in combination with stainless steel components ensures stability for very high strain.

- Extremely small folded size and fast assembly for a high degree of mobility
- Suitable for almost all rigid frame wheelchairs and folding frame models with integrated footrests.
- Can be assembled/disassembled in 30 seconds
- Rolls almost by itself
- No more hot hands downhill!
- Less exertion through better rolling resistance

## Steering head



The unique steering head design delivers razor sharp handling and high steering stability and conceded during cornering and straight-line stability, additionally reinforce two springs this effect!

## Hand-parking brake



Sales price see order sheet.

Color: silver or matt black  
Execution: mechanically

## Frame protection made of plastic



Prevents scratches on the frame of the wheelchair.