

# PELVI.LOC® PL-2DA Active

## PELVI.LOC® PL-2D Thigh Belts

HMV.Nr.: 18.99.99.0905



## Version

**2D-belt for enhanced sports activities with enlarged pads - surface for improved pressure distribution.**

The padded covers give maximum protection of participants and prevent buckle from releasing in use. Hole strap on frame adjustable up to 6 cm.

## Special features

Monoski racer with Martin Braxenthaler:

- **PL-2DA-M** as a chest strap and thigh strap



## Amendment

Tube clamps for attaching the belt to the wheelchair frame available in the following diameters:

## Article numbers

**RS-20:** Tube diameter 20 mm

RS-22: Tube diameter 22 mm

**RS-25:** Tube diameter 25 mm

**RS-30:** Tube diameter 30 mm

**RS-32:** Tube diameter 32 mm



### Scope of application

- Safe driving with 2D pelvic positioning in a wheelchair.
- Can be used for thorax positioning in a standing device.
- Pelvic positioning with pressure distribution and stability in sports wheelchairs.

### Special features

## Article numbers



**PL-2DA-S** for SW 20-32 cm

**PL-2DA-M** for SW 30-45 cm

**PL-2DA-L** for SW 30-45 cm

## Technical Data

	XXXX Bambino	XXS	XS	S	M	L	XL
Seat width (cm)				20 - 32	30 - 40	36 - 46	
Padding size (cm)				19 x 11	23 x 12	28 x 14	
Weight (g)				220	280	280	
Predetermined breaking point (at approx. kg)				60	80	80	

<b>Ratchet size XXXS - XXS</b>			
<b>Ratchet size XS - S</b>	2,4 x 0,2 x 17,0 cm		
<b>Ratchet size M - L - XL</b>	2,4 x 0,2 x 17,0 cm		
<b>Buckle size XXXS</b>			
<b>Buckle size XXS</b>			
<b>Buckle size XS - S</b>	B-24-FLAT-LO-BL		
<b>Buckle size M - L - XL</b>	B-24-ERGO-BL		
<b>Hole-strap</b>	LR-BL-320		