

PELVI.LOC® PL-2DS Basic for moulded Seats

PELVI.LOC® PL-2D Thigh Belts

HMV.Nr.: 18.99.99.0905



Version


2DB-belt suitable for thigh.

The buckles are easily accessible. The ratchet will run in addition to a strap to prevent injuries. Webbing in sizes XS, S, M and L

Amendment

For attachment to the seat shell or on the wheelchair.

Article number
PB-BL-38



Scope of application

- For pelvis positioning in a wheelchair and sports wheelchair.

Article numbers

PL-2DS-XS for width SB 16-24 cm
PL-2DS-S for width SB 20-32 cm
PL-2DS-M for width SB 30-40 cm
PL-2DS-L for width SB 36-46 cm

Technical Data

	XXXS	XXS	XS	S	M	L	XL
Seat width (cm)			16 - 24	20 - 32	30 - 40	36 - 46	
Padding size (cm)			18 x 7	21 x 8	25,5 x 8	26,5 x 8,5	
Weight (g)			220	220	280	280	
Predetermined breaking point (at approx. kg)			60	60	80	80	

Ratchet size XXXS - XXS		Product is latex-free 	 LIKENZIIERTER PARTNER REG. NR. DE 5914274696739-V
Ratchet size XS - S	1,9 x 0,2 x 17,0 cm		
Ratchet size M - L - XL	2,4 x 0,2 x 17,0 cm		
Buckle size XXXS			
Buckle size XXS			
Buckle size XS - S	B-19-ERGO-BL		
Buckle size M - L - XL	B-24-ERGO-BL		