

# PELVI.LOC® PL-3DA Activ

PELVI.LOC® PL-3D Pelvis Belts

HMV.Nr.: 26.99.02.2000



## Version

### 3D-active pelvis retraction belt

3D-belt high tensile forces with enlarged pads - top. Pelvic surface. 3-Dimensional shape gives maximum support of the pelvis without pressure on the spine "anterior superior". The covered buckles help protect from injury and accidental opening during use. Hole strap on frame adjustable up to 6 cm.

## Amendment

Clamping clips for wheelchair frame. Other sizes available on request!

### Article numbers

**RS-20** Tube diameter 20 mm  
**RS-22** Tube diameter 22 mm  
**RS-25** Tube diameter 25 mm  
**RS-30** Tube diameter 30 mm  
**RS-32** Tube diameter 32 mm



## Scope of application



Extensive 3D-belt for high tensile strength (e.g. spasmodic, wheelchair-sports). Covered buckle avoids the unintentional opening of the belt.

### Article numbers

**PL-3DA-XS** for width 16-24 cm  
**PL-3DA-S** for width 20-32 cm  
**PL-3DA-M** for width 30-40 cm  
**PL-3DA-L** for width 36-46 cm

## Technical Data

	XXXS Bambino	XXS	XS	S	M	L	XL
Seat width (cm)			16 - 24	20 - 32	30 - 40	36 - 46	
Padding size (cm)			17,5 x 9	20 x 11	23 x 13	25,5 x 14	
Weight (g)			220	220	280	280	
Predetermined breaking point (at approx. kg)			60	60	80	80	

Ratchet size XXXS - XXS		Product is latex-free 	 LIZENZIIERTER PARTNER REG. NR. DE5914274696739-V
Ratchet size XS - S	1,9 x 0,2 x 17,0 cm		
Ratchet size M - L - XL	2,4 x 0,2 x 17,0 cm		
Buckle size XXXS			
Buckle size XXS			
Buckle size XS - S	B-19-ERGO-BL		
Buckle size M - L - XL	B-24-ERGO-BL		
Hole-strap			