

PELVI.LOC® PL-3DAS Activ for moulded Seat

PELVI.LOC® PL-3D Pelvis Belts

HMV.Nr.: 26.99.02.2000



Version

Pelvis retraction belt for moulded seats

3D-belt high tensile forces with enlarged pads -surface. Nylon strap for passing through the moulded seat. 3-Dimensional shape gives maximum support of the pelvis without pressure on the spine "anterior superior". Formed shape with ultra thin end at the moulded seat fixing point to avoid lateral pressure marks at the pelvis. The covered buckles help protect from injury and accidental opening during use.

Amendment


For moulded seat for better pelvic positioning in the wheelchair.

Article number: PB-BL



Attachment kit for seat buckets.

Article number: BS-24-BL



Scope of application

For custom-made seating systems, moulded seat.

Special features

To improve positioning of the pelvis - retraction of the seat, with high spasticity.

Article numbers

PL-3DAS-XS for width 16-24 cm

PL-3DAS-S for width 20-32 cm

PL-3DAS-M for width 30-40 cm

PL-3DAS-L for width 36-46 cm

PL-3DAS-XL for width 44-52 cm

Technical Data

	XXXS Bambino	XXS	XS	S	M	L	XL
Seat width (cm)			16 - 24	20 - 32	30 - 40	36 - 46	44 - 52
Padding size (cm)			17,5 x 9	20 x 11	23 x 13	28 x 14	29 x 16
Weight (g)			220	220	280	280	280
Predetermined breaking point (at approx. kg)			60	60	80	80	80

Ratchet size XXXS - XXS		Product is latex-free 	
Ratchet size XS - S	1,9 x 0,2 x 17,0 cm		
Ratchet size M - L - XL	2,4 x 0,2 x 17,0 cm		
Ratchet size XL	2,4 x 0,2 x 22,5 cm		
Buckle size XXXS			
Buckle size XXS			
Buckle size XS - S	B-19-ERGO-BL		
Buckle size M - L - XL	B-24-ERGO-BL		
Hole-strap			