

PELVI.LOC® PL-3DB Basic

PELVI.LOC® PL-3D Pelvis Belts

HMV.Nr.: 26.99.02.2000



Version

Pelvis retraction belt


The pelvic supports are 3-dimensional shape for optimal fixation of the pelvis without pressure on the spine "anterior superior ". The ratchet can also be attached with a tab to avoid injury. Hole strap on frame adjustable up to 6 cm.

Amendment

Tube clamps for attaching the belt to the wheelchair frame available in the following diameters.

Article numbers

RS-20 Tube diameter 20 mm
RS-22 Tube diameter 22 mm
RS-25 Tube diameter 25 mm
RS-30 Tube diameter 30 mm
RS-32 Tube diameter 32 mm



Scope of application

Wheelchair, reha-buggy, moulded seat.

Article numbers

PL-3DB-XS for width 16-24 cm
PL-3DB-S for width 20-32 cm
PL-3DB-M for width 30-40 cm
PL-3DB-L for width 36-46 cm
PL-3DB-XL for width 40-48 cm

Technical Data

	XXXS Bambino	XXS	XS	S	M	L	XL
Seat width (cm)			16 - 24	20 - 32	30 - 40	36 - 46	44 - 52
Padding size (cm)			18 x 7	21 x 8	25,5 x 8	26,5 x 8,5	28,5 x 9
Weight (g)			220	220	280	280	280
Predetermined breaking point (at approx. kg)			60	60	80	80	80

Ratchet size XXXS - XXS		Product is latex-free 	 REG. NR. DE5914274696739-V LIZENZIERTER PARTNER
Ratchet size XS - S	1,9 x 0,2 x 17,0 cm		
Ratchet size M - L - XL	2,4 x 0,2 x 17,0 cm		
Buckle size XXXS			
Buckle size XXS			
Buckle size XS - S	B-19-ERGO-BL		
Buckle size M - L - XL	B-24-ERGO-BL		
Hole-strap			