

PELVI.LOC® PL-3DSR Reverso

PELVI.LOC® PL-3D Pelvis Belts

HMV.Nr.: 26.99.02.2000



Version

Formed shape with ultra thin end at the moulded seat fixing point to avoid lateral pressure marks at the pelvis. On the moulded seat remain attached the ratchets. Remove the pads to implement.

Scope of application

For custom made seats, moulded seat.

Special features

The timing belts remain attached to the side of the seat pan. To transfer only remove the pads.

Fine adjustment of the pads to each other can be adjusted via the belt.

Article numbers

PL-3DSR-XXS for SW 12-20 cm

PL-3DSR-XS for SW 16-24 cm

PL-3DSR-S for SW 20-32 cm

PL-3DSR-M for SW 30-40 cm

PL-3DSR-L for SW 36-46 cm

PL-3DSR-XL for SW 44-52 cm

Technical Data

	XXXS Bambino	XXS	XS	S	M	L	XL
Seat width (cm)		12 - 20	16 - 24	20 - 32	30 - 40	36 - 46	44 - 52
Padding size (cm)		15 x 7	18 x 7	21 x 8	25,5 x 8	26,5 x 8,5	28,5 x 9
Weight (g)		220	220	220	280	280	280
Predetermined breaking point (at approx. kg)		60	60	60	80	80	80

Ratchet size XXXS - XXS	1,9 x 0,2 x 40,0 cm	Product is latex-free 	
Ratchet size XS - S	1,9 x 0,2 x 40,0 cm		
Ratchet size M - L - XL	2,4 x 0,2 x 40,0 cm		
Buckle size XXXS			
Buckle size XXS	B-19-ERGO-BL		
Buckle size XS - S	B-19-ERGO-BL		
Buckle size M - L - XL	B-24-ERGO-BL		
Hole-strap XXS - XS - S	V-LR-19-BL-110 and V-LR-19-BL-140		
Hole-strap M - L - XL	V-LR-24-BL-110 and V-LR-24-BL-140		