

PELVI.LOC® PL-3DSRR Reverso with ratchets

PELVI.LOC® PL-3D Pelvis Belts

HMV.Nr.: 26.99.02.2000



Version

Pelvis-retraction-belt with reversed ratchet

Formed shape with ultra thin end at the moulded seat fixing point to avoid lateral pressure marks at the pelvis.

The 3D formed shape is bolted firmly to the timing belt and strap hole.

Buckle belt **R-19/24-BL-400** and pressure-buckle **B-19/24-Star** are included.

Amendment

Hole-strap - set, available in a set of two

- Width: 19/24 mm
- Length: 115 mm

Hole-strap with buckle and clamping clip for mounting onto the frame. Also available in different sizes and diameters.

Example:

- Hole – strap 24 mm wide x 115 mm long and B-24-ERGO and clamping clip RS-25

Article numbers

LR-R-19-BL/B-19-ERGO-TR

LR-R-24-BL/B-24-ERGO-TR



Scope of application

For custom made seats and moulded seat.

Special features

The buckles are mounted on the side of the moulded seat or buggy. The pelvis-retraction-belt with the ratchet slides sideways into the buckles.

To improve positioning of the belt, adjustments can be made via the hole-strap.

Article numbers

PL-3DSRR-XS for width 16-24 cm

PL-3DSRR-S for width 20-32 cm

PL-3DSRR-M for width 30-40 cm

PL-3DSRR-L for width 36-46 cm

Technical Data

	XXXS Bambino	XXS	XS	S	M	L	XL
Seat width (cm)			16 - 24	20 - 32	30 - 40	36 - 46	
Padding size (cm)			18 x 7	21 x 8	25,5 x 8	26,5 x 8,5	
Weight (g)			220	220	280	280	
Predetermined breaking point (at approx. kg)			60	60	80	80	

Ratchet size XXXS - XXS		Product is latex-free 	
Ratchet size XS - S	1,9 x 0,2 x 40,0 cm		
Ratchet size M - L - XL	2,4 x 0,2 x 40,0 cm		
Buckle size XXXS			
Buckle size XXS			
Buckle size XS - S	B-19-STAR-BL		
Buckle size M - L - XL	B-24-STAR-BL		
Hole-strap			