

PELVI.LOC® PL-3DWP Waterproof

PELVI.LOC® PL-3D Pelvis Belts

HMV.Nr.: 26.99.02.2000



Version

Closed-pore special material, impermeable and water repellent! The pelvic pads are 3-dimensionally shaped for optimal fixation of the pelvis without pressure on the anterior superior spine. With easy-to-use stainless steel ratchet buckles, the fit of the strap is adjusted. The perforated belts on the frame can be adjusted up to 8cm in length and thus provide variable space for different situations. To secure the wheelchair user, the pelvic pads are connected in the middle via a toothed belt.



PL-3DWP-WH-L: 36-46 cm, weiß/white



PL-3DWP-BL-L: 36-46 cm, schwarz/black

Amendment

Pipe clamps for attaching the belt to the wheelchair frame available in the following diameters

Article numbers

RS-20 Tube diameter 20 mm
RS-22 Tube diameter 22 mm
RS-25 Tube diameter 25 mm
RS-30 Tube diameter 30 mm
RS-32 Tube diameter 32 mm



Scope of application

For pelvic positioning in the shower wheelchair, bathing chair, active wheelchair, rehab buggy, seat shells and much more.

Article numbers, black

PL-3DWP-BL-XXS for SW 12-20 cm
PL-3DWP-BL-XS for SW 16-24 cm
PL-3DWP-BL-S for SW 20-32 cm
PL-3DWP-BL-M for SW 30-40 cm
PL-3DWP-BL-L for SW 36-46 cm

Article numbers, white

PL-3DWP-WH-XXS for SW 12-20 cm
PL-3DWP-WH-XS for SW 16-24 cm
PL-3DWP-WH-S for SW 20-32 cm
PL-3DWP-WH-M for SW 30-40 cm
PL-3DWP-WH-L for SW 36-46 cm

Special feature

Facet-like thinning of the padding material in the transition from the pad to the fastening strap. Prevents lateral pressure points on the pelvis.

Ratchets, black

V-R-19-BL-ERGO-170
V-R-24-BL-ERGO-170

Ratchets, white

V-R-19-WH-ERGO-170
V-R-24-WH-ERGO-170

Technical Data

	XXXXS Bambino	XXS	XS	S	M	L	XL
Seat width (cm)		12 - 20	16 - 24	20 - 32	30 - 40	36 - 46	44 - 52
Padding size (cm)		14 x 6	18 x 7	21 x 8	25,5 x 8	26,5 x 8,5	28,5 x 9
Weight (g)		220	220	220	280	280	280
Predetermined breaking point (at approx. kg)		40	60	60	80	80	80

Ratchet size XXXS - XXS		1,9 x 0,2 x 17,0 cm	Product is latex-free 	 LIZENZIERTER PARTNER REG. NR. DE 59142746/696739-V
Ratchet size XS - S		1,9 x 0,2 x 17,0 cm		
Ratchet size M - L - XL		2,4 x 0,2 x 17,0 cm		
Buckle size XXXS				
Buckle size XXS				
Buckle size XS - S		B-19-ERGO-BL		
Buckle size M - L - XL		B-24-ERGO-BL		
Hole-strap				