

# PELVI.LOC® PL-CSS Pelvi Contour Seating System

PELVI.LOC® Seating systems (seat trousers)

HMV.Nr.: 26.99.02.2000



## Version

**Pelvic retraction seating system.**

**Colour:** black

The resulting thrust - pressure of the pelvis forward is of the 3D-Contour taken the pads. The padded seat prevents pressure.

## Amendment

Easy and fast assembling the Zahriemens and webbing on the frame with a cable tie, belt loop or pipe clamp.



## Scope of application

For pelvic positioning in a wheelchair or multifunction - Wheelchair.

### Special features

Prevents forward sliding of the basin, especially in reclining backrests and tilt.

Significant reduction in the pressure step by the 3D-Contour of the supports.

Easy handling due to pressure buckle **B-24-Star** star and Timing **R-24-BL-ERGO-400**.

Ratchet and buckle covers, prevents injuries.

### Article numbers

**Pelvi-Contour-Seating-System - Gr. S**

for width 26 cm to 35 cm

**Pelvi-Contour-Seating-System - Gr. M**

for width 36 cm to 39 cm

**Pelvi-Contour-Seating-System - Gr. L**

for width 40 cm to 46 cm

## Technical Data

	XXXS	XXS	XS	S	M	L	XL
Seat width (cm)				26 - 35	36 - 39	40 - 46	
Padding size (cm)				19 x 7,5	29 x 11	36 x 13,5	
Weight (g)				260	300	400	
Predetermined breaking point (at approx. kg)				60	80	80	

Ratchet size XXXS - XXS		Product is latex-free 	 LIZENZIERTER PARTNER REG. NR. DE5914274696739-V
Ratchet size XS - S	2,4 x 0,2 x 40,0 cm		
Ratchet size M - L - XL	2,4 x 0,2 x 40,0 cm		
Buckle size XXXS			
Buckle size XXS			
Buckle size XS - S	B-24-STAR-BL		
Buckle size M - L - XL	B-24-STAR-BL		