

PELVI.LOC® PL-OPS Upper Body Positioning made to measure

PELVI.LOC® / Torso.Loc Upper Body Positioning



HMV.Nr.: 18.99.99.0906

Version

Our upper body positioning system is individually adaptable. According to your measurements, we manufacture the system for safe sitting with upper body stabilization. For the car transport of children, adolescents and adults, the system offers simple and quickly applicable handling, and the upper body positioning can be installed in just a few simple steps.

OPS made to measure

assembled in a few easy steps!

Amendment

Safety buckle protects against unintentional opening of the body belt.



Article number: SF-2

Segufix magnetic closure with 2 magnetic closures incl. 2 security keys to open the system. Only on request and agreement possible!

Article number:

KSS-S-50 3 finger buckle

KSS-50 2 finger buckle

SKS-50-SF Lock with safety function (pin opening)

Scope of application

For the upper body - positioning in a wheelchair, powered wheelchair, car Seat.

Special features

Upper body - positioning - system to your specifications

- Body part with lateral stability
- Seat - belt individually adjustable in length by Velcro fastener
- Wide, padded shoulder straps
- Special – Hook/Loop fasteners with very high tensile force and Life

Article number: PL-OPS-Dimension



Article number: KSS-S-50

To open lock with three fingers - Nylon - Pressure



Article number: KSS-50

To open lock with two fingers - Nylon - Pressure



Article number: SKS-50-SF

Metal lock 50mm wide with safety function (Pin opening)