

PELVI.LOC® PL-OPS made to measure KID

PELVI.LOC® / Torso.Loc Upper Body Positioning

HMV.Nr.: 18.99.99.0906



Version

Our upper body positioning system **KID** is individually adaptable. According to your measurements, we manufacture the system for safe sitting with upper body stabilization. For the car transport of children, adolescents and adults, the system offers simple and quickly applicable handling, and the upper body positioning can be installed in just a few simple steps.

OPS made to measure

assembled in a few easy steps!

Amendment

Various safety buckles for different requirements. To make it more difficult to open the waist belt unintentionally.



Article number:
KSS-S-50 3 finger buckle
KSS-50 2 finger buckle
SKS-30-SF Lock with safety function (pin opening)

Scope of application

For upper body positioning and stabilisation in the car seat.

Special features

PL-OPS-KID consisting of:

- Padded shoulder straps
- Body section with lateral stabilisation
- Seat trousers
- Fastening strap for the backrest of the child seat
- Seat wedge required depending on the age of the child
- Always secure additionally with the 3-point belt in the car!!!

Article number: PL-OPS-KID




Article number: KSS-S-50

To open lock with three fingers - Nylon - Pressure



Article number: KSS-50

To open lock with two fingers - Nylon - Pressure



Article number: SKS-30-SF

Metal lock 30 mm wide with safety function (Pin opening)