

PELVI.LOC® PL-RR Instep Strap

PELVI.LOC® Foot Positioning

HMV.Nr.: 18.99.99.0305

Foot straps

Padded instep strap for securing feet to wheelchair foot plate with 24 mm ratchet and serrated band.



Special features

Delivery always in pairs!
In many buckles and toothed belt variants possible!

Article numbers


- **PL-RR-XS**: Length min. 19 cm **B14-GO-BL** and **R14-GO-110**, Length max. 27 cm
- **PL-RR-S**: Length min. 20 cm **B24-LU-JU** and **R24-BL-170**, Length max. 27 cm
- **PL-RR-M**: Length min. 24 cm **B24-LU-JU** and **R24-BL-195**, Length max. 31 cm
- **PL-RR-L**: Length min. 24 cm **B24-LU-JU** and **R24-BL-225**, Length max. 36 cm
- **PL-RR-XL**: Length min. 29 cm **B24-LU-JU** and **R24-BL-225**, Length max. 51 cm
- **PL-RR-GO-S**: Length min. 20 cm **B19-GO-BL** and **R19-GO-BL-350**, Length max. 27 cm
- **PL-RR-GO-M**: Length min. 24 cm **B19-GO-BL** and **R19-GO-BL-350**, Length max. 31 cm

Included

BS-BL-25-SET

Mounting Bracket Set with screws and nuts.

Always in pairs!



Scope of application

Padded instep strap for securing feet to wheelchair foot plate with 24 mm ratchet and serrated band.
Allows for quick release and securing with one hand.



Version PL-RR

Upholstery size above everything:

- **Size: XS-S** 17,5 cm x 5,0 cm, overall length 21,0 cm
- **Size: M** 17,5 cm x 5 cm, overall length 23,5 cm
- **Size: L** 17,5 cm x 5 cm, overall length 23,5 cm, same cushion size wie Size M



Version PL-RR

Upholstery size above everything:

- **Size: XL** 24,0 cm x 7,0 cm, overall length 35,0 cm



Version PL-RR-GO

Upholstery size above everything:

- **Size: S** 20,0 cm, overall length 27,0 cm
- **Size: M** 24,0 cm, overall length 31,0 cm

Technical data

PL-RR	XS	S	M	L	XL
Padding size (cm)	17,5 x 5	17,5 x 5	25,5 x 8	25,5 x 8	24 x 7
Total length (cm) minimal	19	20	24	24	29
Total length (cm) maximal	27	27	31	36	51

Product is latex free





Technik die verbindet
Technology that connects

Toothed strap	R-14-GO-110	R-24-BL-170	R-24-BL-195	R-24-BL-225	R-24-BL-225
Buckle	B-14-GO-BL	B-24-LU-JU	B-24-LU-JU	B-24-LU-JU	B-24-LU-JU

PL-RR-GO	S	M		
Padding size (cm)	17,5 x 5	25,5 x 8		
Total length (cm) minimal	20	24		
Total length (cm) maximal	27	31		
Toothed strap	R-19-GO-BL-350	R-19-GO-BL-350		
Buckle	B-19-GO-BL	B-19-GO-BL		

Product is latex free